



Viewing Everyday Objects with a Microscope


Materials:

- Microscope
- Everyday objects



Note: *The images provided can be used if a microscope is not available.*

1. Using a microscope, observe daily objects under a microscope.
 - a. Note if an object may be too large, and how a scientist might take a sample.
 - b. Examples include dirt, grass, pencil or eraser shavings, or fingernail clippings.
2. Record observations and make comparisons. Note qualities like if a soft object appears rough, or vice versa.

Glitter



Pencil and shavings



Human fingernail



Eraser and shavings




Penny



Paper



Thread



Eyelash



Computer Screen

